

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Mixture  
 Trade name : Prolystica® Ultra Concentrate HP Neutral Detergent  
 Product code : 1C15

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture : Neutral Detergent

#### 1.3. Details of the supplier of the safety data sheet

STERIS Corporation  
 P. O. Box 147, St. Louis, MO 63166, US  
 Telephone Number for Information: 1-800-548-4873 (Customer Service-Healthcare Products)

#### 1.4. Emergency telephone number

Emergency number : US Emergency Telephone No.1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### GHS-US classification

Oral Tox. 5 H303  
 Derm Tox. 5 H313  
 Skin Irrit. 2 H315  
 Eye Irrit. 2B H320

#### 2.2. Label elements

##### GHS-US labelling

Hazard pictograms (GHS-US) :



GHS07

Signal word (GHS-US) : Warning

Hazard statements (GHS-US) : H303 – May be harmful if swallowed.  
 H313 – May be harmful in contact with skin.  
 H315 – Causes skin irritation.  
 H320 – Causes eye irritation.

Precautionary statements (GHS-US) : P264 – Wash hands thoroughly after handling.  
 P280 – Wear eye protection, protective clothing, protective gloves.  
 P302+P352 – IF ON SKIN: Wash with plenty of soap and water.  
 P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P332+P313 – If skin irritation occurs: Get medical attention.  
 P337+P313 – If eye irritation persists: Get medical attention.  
 P362+P364 – Take off contaminated clothing and wash before reuse.

#### 2.3. Other hazards

No additional information available.

#### 2.4. Unknown acute toxicity (GHS-US)

No data available.

### SECTION 3: Composition/information on ingredients

#### 3.1. Substance

Not applicable.

# Prolystica<sup>®</sup> Ultra Concentrate HP Neutral Detergent

## Safety Data Sheet

according to the federal final rule of hazard communication revised on 2012 (HazCom 2012)

### 3.2. Mixture

Name	Product identifier	%	GHS-US classification
Tolyltriazole, sodium salt	(CAS No) 64665-57-2	1 - 5	Acute Tox. 4 (Oral), H302 Skin Corr. 1B, H314 Eye Dam. 1, H318 Aquatic Chronic 2, H411
Disodium Ethylhexyliminodipropionate	(CAS No) 68610-44-6	1 - 5	Not classified
Octanaminium, N-(Carboxymethyl) – N, N-Dimethyl-Inner Salt	(CAS No) 27593-14-2	7 - 13	Not classified

Full text of H-phrases: see Section 16.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general	: Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove patient to fresh air and keep at rest in a position comfortable for breathing. If not breathing, give artificial respiration. Get medical attention.
First-aid measures after skin contact	: Immediately flush skin with plenty of water for at least 15 minutes. Remove/Take off immediately all contaminated clothing. Obtain medical attention.
First-aid measures after eye contact	: In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist. Remove contact lenses, if present and easy to do. Continue rinsing. Seek medical attention immediately.
First-aid measures after ingestion	: If swallowed, rinse mouth with water (only if the person is conscious). Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting. Give water to drink if victim completely conscious/alert.

### 4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries after inhalation	: The inhalation of airborne droplets or aerosols causes irritation of the respiratory tract.
Symptoms/injuries after skin contact	: Causes skin irritation.
Symptoms/injuries after eye contact	: Causes eye irritation.
Symptoms/injuries after ingestion	: May cause gastrointestinal tract irritation.

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media	: Use fire-extinguishing media appropriate for surrounding materials. Foam, carbon dioxide, dry chemical. Water spray.
------------------------------	--

### 5.2. Special hazards arising from the substance or mixture

No additional information available.

### 5.3. Advice for firefighters

Firefighting instructions	: Prevent fire-fighting water from entering environment. Cool closed containers exposed to fire with water spray. Do not get water inside containers.
Protective equipment for firefighters	: Use self-contained breathing apparatus. Do not enter fire area without proper protective equipment, including respiratory protection.
Other information	: Hazardous decomposition products may be released during prolonged heating like smokes, carbon monoxide and dioxide.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

General measures	: Avoid inhalation of vapor and spray mist. Avoid contact with skin, eyes and clothes. Use personal protective equipment as required. Stop leak if safe to do so.
------------------	---

#### 6.1.1. For non-emergency personnel

Emergency procedures	: Evacuate unnecessary personnel.
----------------------	-----------------------------------

#### 6.1.2. For emergency responders

Protective equipment	: Equip cleanup crew with proper protection.
Emergency procedures	: Ensure adequate ventilation.

# Prolystica® Ultra Concentrate HP Neutral Detergent

## Safety Data Sheet

according to the federal final rule of hazard communication revised on 2012 (HazCom 2012)

### 6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters. Avoid release to the environment.

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Spills may be picked up with a mop and followed by a water rinse. Wash contaminated areas with large quantities of water to a sanitary sewer, if in accordance with local, state or national legislation. Large spills: Contain any spills with dikes or absorbents to prevent migration and entry into sewers or streams. Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible. Collect in closed containers for disposal. Ensure all national/local regulations are observed.

### 6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Read label before use. Avoid contact with skin, eyes and clothing. Avoid breathing mist or vapor . Keep container closed when not in use. Work in a well-ventilated area.

Hygiene measures : Wash hands thoroughly after handling. Take care for general good hygiene and housekeeping. Do not eat, drink or smoke when using this product. Wash contaminated clothing before reuse.

### 7.2. Conditions for safe storage, including any incompatibilities

Technical measures : A washing facility/water for eye and skin cleaning purposes should be present. Provide local exhaust or general room ventilation.

Storage conditions : Keep out of reach of children. Keep only in the original container in a cool, well ventilated place.

Storage temperature : When dispensing, store in a cool place below 32°C (90°F) to maintain product integrity

### 7.3. Specific end use(s)

No additional information available.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available.

### 8.2. Exposure controls

Appropriate engineering controls : Ensure adequate ventilation. Control airborne concentrations below the exposure limits. Emergency eye wash fountains and safety showers should be available in the immediate vicinity of any potential exposure.

Personal protective equipment : Avoid all unnecessary exposure. Personal protective equipment should be selected based upon the conditions under which this product is handled or used. Protective clothing. Gloves. Protective goggles.



Hand protection : Wear rubber gloves.

Eye protection : Wear chemical goggles or safety glasses.

Respiratory protection : In case of insufficient ventilation, wear suitable respiratory equipment. Where exposure through inhalation may occur from use, respiratory protection equipment is recommended.

Other information : Do not eat, drink or smoke during use.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state : Liquid

Color : Yellow to amber

Odor : Mild chemical odor

Odor threshold : No data available

pH : 7.7

Relative evaporation rate (butyl acetate=1) : No data available

Melting point : No data available

Freezing point : No data available

# Prolystica<sup>®</sup> Ultra Concentrate HP Neutral Detergent

## Safety Data Sheet

according to the federal final rule of hazard communication revised on 2012 (HazCom 2012)

Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Density	: ca. 1.12 g/ml Specific Gravity
Solubility	: Water: Completely soluble
Log Pow	: No data available
Log Kow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available.

### 9.2. Other information

No additional information available.

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No additional information available.

### 10.2. Chemical stability

Stable under normal conditions of use.

### 10.3. Possibility of hazardous reactions

Hazardous polymerization will not occur.

### 10.4. Conditions to avoid

No additional information available.

### 10.5. Incompatible materials

No additional information available.

### 10.6. Hazardous decomposition products

No additional information available.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified Based on available data, the classification criteria are not met.
----------------	---

#### Tolytriazole, sodium salt (64665-57-2)

ATE CLP (oral)	500.000 mg/kg bodyweight
Skin corrosion/irritation	: Causes skin irritation pH: 7.7
Serious eye damage/irritation	: Irritating to eyes pH: 7.7
Respiratory or skin sensitisation	: Not classified Based on available data, the classification criteria are not met
Germ cell mutagenicity	: Not classified Based on available data, the classification criteria are not met
Carcinogenicity	: Not classified Based on available data, the classification criteria are not met

# Prolystica<sup>®</sup> Ultra Concentrate HP Neutral Detergent

## Safety Data Sheet

according to the federal final rule of hazard communication revised on 2012 (HazCom 2012)

Reproductive toxicity	: Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (single exposure)	: Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (repeated exposure)	: Not classified Based on available data, the classification criteria are not met
Aspiration hazard	: Not classified Based on available data, the classification criteria are not met.

### SECTION 12: Ecological information

#### 12.1. Toxicity

Not applicable.

#### 12.2. Persistence and degradability

No additional information available.

#### 12.3. Bioaccumulative potential

##### Prolystica<sup>®</sup> Ultra Concentrate HP Neutral Detergent

Bioaccumulative potential	Not established.
---------------------------	------------------

#### 12.4. Mobility in soil

No additional information available.

#### 12.5. Other adverse effects

Other information : Avoid release to the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

Waste disposal recommendations	: Dispose in a safe manner in accordance with local/national regulations. Do not re-use empty containers.
Additional information	: Never return unused material to original container. Small spills may be flushed to a sanitary sewer with copious amounts of water, if in accordance with local, state or national legislation. Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials	: Avoid release to the environment.

### SECTION 14: Transport information

In accordance with DOT  
Not regulated for transport.

#### Additional information

Other information : No supplementary information available.

#### ADR

#### Transport by sea

No additional information available.

#### Air transport

No additional information available.

### SECTION 15: Regulatory information

#### 15.1. US Federal regulations

##### Prolystica<sup>®</sup> Ultra Concentrate HP Neutral Detergent

RQ (Reportable quantity, section 304 of EPA's List of Lists) :	1250000 lb
--	------------

##### 1-Octanamine, N,N-dimethyl- (7378-99-6)

Listed on the United States TSCA (Toxic Substances Control Act) inventory

##### Tolyltriazole, sodium salt (64665-57-2)

Listed on the United States TSCA (Toxic Substances Control Act) inventory

# Prolystica<sup>®</sup> Ultra Concentrate HP Neutral Detergent

## Safety Data Sheet

according to the federal final rule of hazard communication revised on 2012 (HazCom 2012)

### 15.2. International regulations

Not applicable.

### 15.3. US State regulations

Not applicable.

## SECTION 16: Other information

Revision Date : 03/26/2015

Full text of H-phrases:

Acute Tox. 4 (Oral)	Acute toxicity (Oral), Category 4
Acute Tox. 5 (Derm)	Acute toxicity (Derm), Category 5
Acute Tox. 5 (Oral)	Acute toxicity (Oral), Category 5
Aquatic Chronic 2	Hazardous to the aquatic environment — Chronic Hazard, Category 2
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Eye Irrit. 2B	Serious eye damage/eye irritation, Category 2B
Skin Corr. 1B	Skin corrosion/irritation Category 1B
Skin Irrit. 2	Skin corrosion/irritation Category 2
H302	Harmful if swallowed
H303	May be harmful if swallowed
H313	May be harmful in contact with skin
H314	Causes severe skin burns and eye damage
H315	Causes skin irritation
H318	Causes serious eye damage
H320	Causes eye irritation
H411	Toxic to aquatic life with long lasting effects

NFPA health hazard

: 1 - Exposure could cause irritation but only minor residual injury even if no treatment is given.

NFPA fire hazard

: 0 - Materials that will not burn.

NFPA reactivity

: 0 - Normally stable, even under fire exposure conditions, and are not reactive with water.



SDS US (GHS HazCom 2012)

*The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.*