

**1. Identification of the Substance and Company**

**Chemdi® Indicator Strips**

**Product No.**

**SDS No. S111UK**

**Prepared by: M. Ebers**

[asksteris\\_msds@steris.com](mailto:asksteris_msds@steris.com)

**Date Created: 12/01/2001**

**Date Revised: 01/23/2018**

**Date Reviewed: NA**

**NFPA 704 HAZARD RATING:**

**HEALTH: 0**

**FIRE: 0**

**REACTIVITY: 0**

STERIS Corporation, P. O. Box 147, St. Louis, MO 63166, US  
Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)  
Telephone Number for Information: 1-800-548-4873 (Customer Service-Healthcare Products)

STERIS Ireland Limited, IDA Business and Technology Park, Tullamore  
County Offaly, R35 X865, Ireland.  
Emergency Phone No: +44 (0) 1895 622 639  
Product/Technical Information Phone No: +44 (0) 116 276

<b>PRODUCT NAME:</b>	<b>PRODUCT NUMBER:</b>
Steam – 8½" – Standard Pack	NB001
Steam – 8½" – Bulk Pack	NB108
Steam – 4" – Standard Pack	NB034
EO- 8½" – Standard Pack	NB002
EO – 8½" – Bulk Pack	NB110
EO- 4" – Standard Pack	NB035
Dual Indicator Strips (EO and Steam)	SC705
<b>SHIPPING INFORMATION:</b>	<b>NOT RESTRICTED</b>

Chemdi® Indicator Strips are NOT hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200.

Chemdi® Indicator Strips are designed to monitor achievement of sterilization cycle parameters for steam and EO sterilization. NOT to be used as a biological indicator to assure sterility.