

1. Identification of the Substance and Company

Temp-A-Sure™ Indicator Strips

Product No.

SDS No. S114UK

Prepared by: M. Ebers

asksteris_msd@steris.com

Date Created: 12/01/2001

Date Revised: 02/08/2018

NFPA 704 HAZARD RATING:

HEALTH: 0

FIRE: 0

REACTIVITY: 0

Date Reviewed: NA

STERIS Corporation, P. O. Box 147, St. Louis, MO 63166, US
Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)
Telephone Number for Information: 1-800-548-4873 (Customer Service-Healthcare Products)

STERIS Ireland Limited, IDA Business and Technology Park, Tullamore
County Offaly, R35 X865, Ireland.
Emergency Phone No: +44 (0) 1895 622 639
Product/Technical Information Phone No: +44 (0) 116 276 8636

PRODUCT NAME:	PRODUCT NUMBER:
Temp-A-Sure Indicator Strips	NB225
SHIPPING INFORMATION:	NOT RESTRICTED

Temp-A-Sure Indicator Strips are NOT hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200.

Temp-A-Sure Indicator Strips measure water temperature in the washer. Four windows 160°F (71°C), 170°F (77°C), 180°F (83°C) and 190°F (88°C) indicate process temperature is achieved.