

**1. Identification of the Substance and Company**

**VERIFY™ Bowie Dick Test Pack**

Product No EQC009 and EQC010

SDS No. B102UK

Prepared by: Richard Bancroft

[asksteris\\_msds@steris.com](mailto:asksteris_msds@steris.com)

Date Created: 12/10/2004

Date Revised: 02/12/2018

**NFPA 704 HAZARD RATING:**

HEALTH: 0

FIRE: 0

REACTIVITY: 0

Date Reviewed: NA

STERIS Corporation, P.O. Box 147, St. Louis, Mo 63166, US  
Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)  
Telephone Number for information: 1-800-548-4873 (Customer Service – Healthcare Products)

STERIS Ireland Limited, IDA Business and Technology Park, Tullamore  
County Offaly, R35 X865, Ireland.  
Emergency Phone No: +44 (0) 1895 622 639  
Product/Technical Information Phone No: +44 (0) 116 276

VERIFY™ Bowie Dick Test Packs are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200 and the Hazardous Products Regulations; and are free of lead and other heavy metals.

VERIFY™ Bowie Dick Test Packs do not contain any components derived from animal origin, do not come into contact with materials of animal origin and therefore comply with the "Note for guidance of minimizing the risk of transmitting animal spongiform encephalopathy agents via human Veterinary medicinal products (EMEA/410/01 Rev. 2-October 2003)."

VERIFY™ Bowie Dick Test Packs contain natural dry rubber.

The VERIFY™ Bowie Dick Test Pack verifies that the prevacuum sterilizer effectively removes air and provides efficient steam penetration. The advanced thermochromic ink formulation can detect problems with steam quality.

*The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.*