



1. Identification of the Substance and Company
VERIFY™ Valueline STEAM Indicator Strips

Product No. PCC009
SDS No. V100US

NFPA 704 HAZARD RATING:

HEALTH: 0
FIRE: 0
REACTIVITY: 0

Prepared by: Richard Bancroft
asksteris_msds@steris.com

Date Created: 07/26/2006

Date Revised: 02/19/2018

Date Reviewed: N/A

STERIS Corporation, P.O. Box 147, St. Louis, Mo 63166, US
Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)
Telephone Number for information: 1-800-548-4873 (Customer Service – Healthcare Products)

STERIS Ireland Limited, IDA Business and Technology Park, Tullamore
County Offaly, R35 X865, Ireland.
Emergency Phone No: +44 (0) 1895 622 639
Product/Technical Information Phone No: +44 (0) 116 276 8636

Temperature:	Product Number
250-275°F (121-135°C)	PCC009
SHIPPING INFORMATION:	NOT RESTRICTED

VERIFY™ Valueline **STEAM** Indicator Strips are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200 and the Controlled Products Regulations, and are manufactured in the absence of lead and other heavy metals.

VERIFY™ Valueline **STEAM** Indicator Strips are manufactured in the absence of latex and natural dry rubber products.

VERIFY™ Valueline **STEAM** Indicator Strips are Type 1 process indicators as defined by ANSI / AAMI / ISO 11140-1 designed to confirm that a device has been exposed to a steam sterilization cycle.

VERIFY™ Valueline **STEAM** Indicator Strips do not contain any components derived from animal origin, do not come into contact with materials of animal origin and therefore comply with the "Note for guidance of minimizing the risk of transmitting animal spongiform encephalopathy agents via human Veterinary medicinal products (EMEA/410/01 Rev. 2-October 2003)."

The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.