



**STERIS**

**Spor dex® VH2O2**

**BIOLOGICAL INDICATOR MEDIA**

Safety Data Sheet

Date of issue: March 7, 2014

Version: 1.0

**1. Identification of the Substance and Company**

**Spor dex® VH2O2 Biological Indicator Media**

**Product No.** NA117

**SDS No.** S291

**Prepared by:** M. Ebers

[asksteris\\_msds@steris.com](mailto:asksteris_msds@steris.com)

**Date Created:** 03/07/2014

**Date Revised:** 02/19/2018

**Date Reviewed:** N/A

**NFPA 704 HAZARD RATING:**

**HEALTH:** 0

**FIRE:** 0

**REACTIVITY:** 0

STERIS Corporation, P.O. Box 147, St. Louis, Mo 63166, US  
Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)  
Telephone Number for information: 1-800-548-4873 (Customer Service – Healthcare Products)

STERIS Ireland Limited, IDA Business and Technology Park, Tullamore  
County Offaly, R35 X865, Ireland.  
Emergency Phone No: +44 (0) 1895 622 639  
Product/Technical Information Phone No: +44 (0) 116 276 8636

<b>PRODUCT NAME:</b>	
<b>Spor dex® VH2O2 Biological Indicator Media</b>	
<b>SHIPPING INFORMATION:</b>	<b>NOT RESTRICTED</b>

Spor dex® VH2O2 Biological Indicator Media is NOT considered to be a hazardous product according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200 and Regulation (EC) No. 453/2010.

Spor dex® VH2O2 Biological Indicator Media is composed of Trypticase Soy Broth and a pH indicator. Appropriate gloves and safety glasses should be worn when removing this product from the sterilizer or cleaning spilled media/broken glass.

Store at room temperature. Autoclave at 121°C for 30 minutes prior to disposal, or incinerate.

*The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.*