

**1. Identification of the Substance and Company**
**Spordex® Single-Species Spore Strips and Discs**
**Product No.**
**SDS No. S106UK**
**NFPA 704 HAZARD RATING:**
**HEALTH: 0**
**FIRE: 0**
**REACTIVITY: 0**
**Prepared by: M. Ebers**
[asksteris\\_msds@steris.com](mailto:asksteris_msds@steris.com)
**Date Created: 12/01/2001**
**Date Revised: 02/08/2018**
**Date Reviewed: N/A**

STERIS Corporation, P. O. Box 147, St. Louis, MO 63166, US  
 Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)  
 Telephone Number for Information: 1-800-548-4873 (Customer Service-Healthcare Products)

STERIS Ireland Limited, IDA Business and Technology Park, Tullamore  
 County Offaly, R35 X865, Ireland.

Emergency Phone No: +44 (0) 1895 622 639

Product/Technical Information Phone No: +44 (0) 116 276 8636

<b>PRODUCT NAME:</b>	<b>PRODUCT NUMBER:</b>
<b>For EO and Dry Heat:</b>	
<i>10<sup>8</sup> B. atrophaeus</i> - strips	NA034
<i>10<sup>7</sup> B. atrophaeus</i> - strips	NA033
<i>10<sup>6</sup> B. atrophaeus</i> - strips	NA005
<i>10<sup>6</sup> B. atrophaeus</i> - strips	NA004
<i>10<sup>5</sup> B. atrophaeus</i> - strips	NA032
<i>10<sup>4</sup> B. atrophaeus</i> - strips	NA031
<i>10<sup>3</sup> B. atrophaeus</i> - strips	NA112
<i>10<sup>8</sup> B. atrophaeus</i> - discs	NA082
<i>10<sup>7</sup> B. atrophaeus</i> - discs	NA105
<i>10<sup>6</sup> B. atrophaeus</i> - discs	NA075
<i>10<sup>6</sup> B. atrophaeus</i> - discs	NA130
<i>10<sup>3</sup> B. atrophaeus</i> - discs	NA120
<b>For Steam:</b>	
<i>10<sup>6</sup> G. stearothermophilus</i> - strips (D-value = 1.5-2.5)	NA037
<i>10<sup>6</sup> G. stearothermophilus</i> - strips (D-value = 1.5-2.5)	NA039
<i>10<sup>6</sup> G. stearothermophilus</i> - strips (D-value = 2.0-2.5)	NA141
<i>10<sup>6</sup> G. stearothermophilus</i> - strips (D-value = 2.0-2.5)	NA139
<i>10<sup>6</sup> G. stearothermophilus</i> - strips (D-value = 1.5-1.9)	NA142
<i>10<sup>6</sup> G. stearothermophilus</i> - strips (D-value = 1.5-1.9)	NA143
<i>10<sup>5</sup> G. stearothermophilus</i> - strips	NA036
<i>10<sup>4</sup> G. stearothermophilus</i> - strips	NA007
<i>10<sup>4</sup> G. stearothermophilus</i> - strips	NA006
<i>10<sup>3</sup> G. stearothermophilus</i> - strips	NA035
<i>10<sup>6</sup> G. stearothermophilus</i> - discs	NA111



<b>SHIPPING INFORMATION:</b>
------------------------------

<b>NOT RESTRICTED</b>
-----------------------

The Spordex® products are NOT hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200.

Spordex® strips and discs are inoculated with spores of *Geobacillus stearothermophilus* for steam sterilization, or *Bacillus atrophaeus* for EO and dry heat use.

The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.