

**1. Identification of the Substance and Company**  
**VERIFY<sup>™</sup> SixCess<sup>™</sup> Steam Indicator Strips**

 Product No. See Chart Below  
 SDS No. B111US
**NFPA 704 HAZARD RATING:**
**HEALTH:** 0  
**FIRE:** 0  
**REACTIVITY:** 0

 Prepared by: M. Ebers  
[asksteris\\_msds@steris.com](mailto:asksteris_msds@steris.com)  
 Date of Issue: 06/04/2015

Revision Date: 02/16/2018

 STERIS Corporation, P.O. Box 147, St. Louis, Mo 63166, US  
 Emergency Telephone No. 1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)  
 Telephone Number for information: 1-800-548-4873 (Customer Service – Healthcare Products)

 STERIS Ireland Limited, IDA Business and Technology Park, Tullamore  
 County Offaly, R35 X865, Ireland.  
 Emergency Phone No: +44 (0) 1895 622 639  
 Product/Technical Information Phone No: +44 (0) 116 276 8636

Product Name:	Product Number
VERIFY <sup>™</sup> SixCess <sup>™</sup> 250F 30 min Short Steam Indicator	PCC005
VERIFY <sup>™</sup> SixCess <sup>™</sup> 270F 4 min Short Steam Indicator	PCC006
VERIFY <sup>™</sup> SixCess <sup>™</sup> 250F 30min Long Steam Indicator	PCC007
VERIFY <sup>™</sup> SixCess <sup>™</sup> 270F 4 min Long Steam Indicator	PCC008
VERIFY <sup>™</sup> SixCess <sup>™</sup> 270F 15 min Steam Indicator	PCC031
VERIFY <sup>™</sup> SixCess <sup>™</sup> 275F 10 min Long Steam Indicator	PCC033
VERIFY <sup>™</sup> SixCess <sup>™</sup> 275F 3 min Long Steam Indicator	PCC034
VERIFY <sup>™</sup> SixCess <sup>™</sup> 275F 10 min Short Steam Indicator	PCC037
VERIFY <sup>™</sup> SixCess <sup>™</sup> 275F 3 min Short Steam Indicator	PCC038
VERIFY <sup>™</sup> SixCess <sup>™</sup> 270F Flash Indicator	LCC310
VERIFY <sup>™</sup> SixCess <sup>™</sup> 275F Flash indicator	LCC005
<b>SHIPPING INFORMATION:</b>	<b>NOT RESTRICTED</b>

VERIFY<sup>™</sup> SixCess<sup>™</sup> steam chemical indicators are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200 and the Hazardous Products Regulations; and are free of lead and other heavy metals.

VERIFY<sup>™</sup> SixCess<sup>™</sup> steam chemical indicators do not contain any components derived from animal origin, do not come into contact with materials of animal origin and therefore comply with the "Note for guidance of minimizing the risk of transmitting animal spongiform encephalopathy agents via human Veterinary medicinal products (EMA/410/01 Rev. 2-October 2003)."

The VERIFY<sup>™</sup> SixCess<sup>™</sup> series of steam chemical indicators are free of natural rubber latex and dry natural rubber products.

VERIFY<sup>™</sup> SixCess<sup>™</sup> Steam Indicators are cycle verification indicators which are designed to react to all critical variables for specific sterilization cycles. These cycle-specific indicators have a more defined endpoint (point of color change) than Class 5 (integrating indicators. Each strip monitors all critical parameters of steam sterilization (exposure time, temperature, and steam quality).

*The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.*