

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Mixture  
Trade name : PRE-KLENZ™ Soak Shield  
Product code : 2D66

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture : Professional use  
Use of the substance/mixture : Presoak

#### 1.3. Details of the supplier of the safety data sheet

STERIS Corporation  
P. O. Box 147, St. Louis, MO 63166, US  
Telephone Number for Information: 1-800-548-4873 (Customer Service-Healthcare Products)

Supplier:

Device Technologies Australia Pty Ltd  
1 Garigal Road  
Belrose NSW 2085, Australia  
Telephone: 1-800-429-551  
Fax: 612-9975-5711

Device Technologies New Zealand Limited  
47 Arrenway Drive, Albany, Auckland, 0632  
New Zealand  
Tel: 0508 338 423, Fax: 649 9913 2009.

#### 1.4. Emergency telephone number

Emergency number : US Emergency Telephone No.1-314-535-1395 (STERIS); 1-800-424-9300 (CHEMTREC)  
1 800 429 551 (24 hours) Australia  
0508 338 423 (New Zealand)

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

Not classified as hazardous per GHS-AU and ADG

#### 2.2. Label elements

##### GHS-AU labelling

Not classified

#### 2.3. Other hazards

No additional information available

#### 2.4. Unknown acute toxicity (GHS-AU)

No data available

### SECTION 3: Composition/information on ingredients

#### 3.1. Substance

Not applicable

#### 3.2. Mixture

Not applicable

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).  
First-aid measures after inhalation : Remove patient to fresh air and keep at rest in a position comfortable for breathing. If not breathing, give artificial respiration. Seek medical attention immediately.  
First-aid measures after skin contact : Remove/Take off immediately all contaminated clothing. Immediately flush skin with plenty of water for at least 15 minutes. If skin irritation persists, seek medical attention.

# PRE-KLENZ™ Soak Shield

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

- First-aid measures after eye contact : Immediately flush eyes thoroughly with water for at least 15 minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Seek medical attention immediately.
- First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Give water to drink if victim completely conscious/alert. Immediately call a POISON CENTER or doctor/physician.

### 4.2. Most important symptoms and effects, both acute and delayed

- Symptoms/injuries after inhalation : May cause an allergic skin reaction. Dust of the product, if present, may cause respiratory irritation after an excessive inhalation exposure. May cause allergy or asthma symptoms or breathing difficulties if inhaled. Inhalation may cause irritation, cough, short breathing.
- Symptoms/injuries after skin contact : Repeated or prolonged skin contact may cause irritation.
- Symptoms/injuries after eye contact : May cause irritation.
- Symptoms/injuries after ingestion : Can occur: Gastrointestinal disturbance.

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

- Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire.

### 5.2. Special hazards arising from the substance or mixture

No additional information available

### 5.3. Advice for firefighters

- Firefighting instructions : Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.
- Protective equipment for firefighters : Do not enter fire area without proper protective equipment, including respiratory protection.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

- General measures : Stop leak if safe to do so. Avoid contact with skin, eyes and clothing. Avoid breathing dust, mist or spray. Spilled material may present a slipping hazard. Ensure adequate air ventilation. Work in a well-ventilated area.

#### 6.1.1. For non-emergency personnel

- Emergency procedures : Evacuate unnecessary personnel.

#### 6.1.2. For emergency responders

- Protective equipment : Equip cleanup crew with proper protection.
- Emergency procedures : Ventilate area.

### 6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

### 6.3. Methods and material for containment and cleaning up

- Methods for cleaning up : Spills may be picked up with a mop and followed by a water rinse. Small spills may be flushed to a sanitary sewer with copious amounts of water. Ensure all national/local regulations are observed.

### 6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

- Precautions for safe handling : Avoid all eye and skin contact and do not breathe vapour and mist. Provide good ventilation in process area to prevent formation of vapour. For further information refer to section 8 : Exposure-controls/personal protection.
- Hygiene measures : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Separate working clothes from town clothes. Launder separately. Handle in accordance with good industrial hygiene and safety practices.

# PRE-KLENZ™ Soak Shield

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

### 7.2. Conditions for safe storage, including any incompatibilities

- Technical measures : A washing facility/water for eye and skin cleaning purposes should be present. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of mists and/or vapours below the recommended exposure limits.
- Storage conditions : Keep out of reach of children. Keep only in original container. Keep container closed when not in use. Store in a dry, cool and well-ventilated place. Keep the packing dry and well sealed to prevent contamination and absorption of humidity. Avoid high temperatures.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

- Appropriate engineering controls : Emergency eye wash fountains and safety showers should be available in the event of any potential exposure. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of mists and/or vapors below the recommended exposure limits.
- Personal protective equipment : Avoid all unnecessary exposure. Personal protective equipment should be selected based upon the conditions under which this product is handled or used. Protective clothing. Gloves.



- Hand protection : Wear rubber gloves.
- Eye protection : Wear safety glasses.
- Skin and body protection : Wear suitable protective clothing.
- Respiratory protection : Fine dispersion/spraying/misting: In case of insufficient ventilation, wear suitable respiratory equipment.
- Other information : Do not eat, drink or smoke during use.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

- Physical state : Liquid
- Colour : Amber.
- Odour : Spicy floral fragrance
- Odour threshold : No data available
- pH : 7.5 - 8 Approximately
- Relative evaporation rate (butyl acetate=1) : No data available
- Melting point : No data available
- Freezing point : No data available
- Boiling point : No data available
- Flash point : No data available
- Auto-ignition temperature : No data available
- Decomposition temperature : No data available
- Flammability (solid, gas) : No data available
- Vapour pressure : No data available
- Relative vapour density at 20 °C : No data available
- Relative density : No data available
- Density : 1.000 g/ml Specific Gravity
- Solubility : Water: completely soluble
- Log Pow : No data available
- Log Kow : No data available

# PRE-KLENZ™ Soak Shield

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No additional information available.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

Hazardous polymerization will not occur.

### 10.4. Conditions to avoid

No additional information available.

### 10.5. Incompatible materials

Strong oxidizers.

### 10.6. Hazardous decomposition products

Thermal decomposition generates : Fume. Carbon monoxide. Carbon dioxide. Aldehydes.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified (Based on available data, the classification criteria are not met)
Skin corrosion/irritation	: Not classified Based on available data, the classification criteria are not met. pH: 7.5 - 8 Approximately
Serious eye damage/irritation	: Not classified Based on available data, the classification criteria are not met. pH: 7.5 - 8 Approximately
Respiratory or skin sensitisation	: May cause an allergic skin reaction.
Germ cell mutagenicity	: Not classified Based on available data, the classification criteria are not met.
Carcinogenicity	: Not classified Based on available data, the classification criteria are not met.
Reproductive toxicity	: Not classified Based on available data, the classification criteria are not met.
Specific target organ toxicity (single exposure)	: Not classified Based on available data, the classification criteria are not met.
Specific target organ toxicity (repeated exposure)	: Not classified Based on available data, the classification criteria are not met.
Aspiration hazard	: Not classified Based on available data, the classification criteria are not met.
Potential Adverse human health effects and symptoms	: Not classified Based on available data, the classification criteria are not met.

# PRE-KLENZ™ Soak Shield

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

- Symptoms/injuries after inhalation : May cause an allergic skin reaction. Dust of the product, if present, may cause respiratory irritation after an excessive inhalation exposure. May cause allergy or asthma symptoms or breathing difficulties if inhaled. Inhalation may cause irritation, cough, short breathing.
- Symptoms/injuries after skin contact : Repeated or prolonged skin contact may cause irritation.
- Symptoms/injuries after eye contact : May cause moderate irritation.
- Symptoms/injuries after ingestion : Can occur: Gastrointestinal disturbance.

## SECTION 12: Ecological information

### 12.1. Toxicity

- Ecology - general : Readily biodegradable. (OECD 301B method). The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No. 648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer.

### 12.2. Persistence and degradability

#### PRE-KLENZ Soak Shield

Persistence and degradability	Not established.
-------------------------------	------------------

### 12.3. Bioaccumulative potential

#### PRE-KLENZ Soak Shield

Bioaccumulative potential	Not established.
---------------------------	------------------

### 12.4. Mobility in soil

No additional information available

### 12.5. Other adverse effects

- Other information : Avoid release to the environment.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

- Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Do not allow to enter into surface water or drains. Ensure all national/local regulations are observed. Empty containers should be rinsed and discarded in a waste container or offered for recycling if possible.
- Ecology - waste materials : Avoid release to the environment.

## SECTION 14: Transport information

- In accordance with DOT  
Not regulated for transport  
Not classified as hazardous per GHS-AU and ADG

### Additional information

- Other information : No supplementary information available.

### ADR

- Transport document description : No additional information available

### Transport by sea

No additional information available

### Air transport

No additional information available

## SECTION 15: Regulatory information

### 15.1. US Federal regulations

Not applicable

#### 15.1.1. Australia

AICS Listed or Exempt, Not classified as hazardous per GHS-AU or ADG

# PRE-KLENZ™ Soak Shield

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

### 15.2. International regulations

Not applicable

### 15.3. US State regulations

Not applicable.

## SECTION 16: Other information


Revision Date : 06/13/2018

Other information : None.

NFPA health hazard : 1 - Exposure could cause irritation but only minor residual injury even if no treatment is given.

NFPA fire hazard : 0 - Materials that will not burn.

NFPA reactivity : 0 - Normally stable, even under fire exposure conditions, and are not reactive with water.



SDS US (GHS HazCom 2012)

*The information on this sheet is not a specification and does not guarantee specific properties. The information is intended to provide general knowledge as to health and safety based upon our knowledge of the handling, storage and use of the product. It is not applicable to unusual or non-standard uses of the product or where instruction or recommendations are not followed.*